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QUICK-GUIDE

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GM Counter w/ built-in tube 5135.45

QUICK GUIDE

Power On: Press **ON** to power on the instrument. The display briefly shows a version number and the GM voltage. Immediately after the counter has been turned on, the counting period is 1 s – **CPS** is indicated on the display.

Other Views: From **CPS** press **▲** once to select **μSv/h** relative to Cs-137.

Press **▲** once more to select the mean value of counts per second – **Mean CPS**.

Other counting periods: Press **▲** or **▼** until **Pulse** is highlighted in the display. The instrument now counts for a selectable period, starting with 10 s.

Press **HOLD ▲ HOLD** and set the measurement period with **▲** and **▼**.

The unit is seconds. The symbol "----" means indefinite counting.

Select with **CAL** / cancel with **HOLD**.

Setting the GM voltage: Press **CAL** and adjust the voltage with **▲** and **▼**.

Select with **CAL** / cancel with **HOLD**.

Counting rate (**CPS**) is shown in the secondary display during the adjustment.

This setting returns to its default value when the instrument is turned off.

- 8 fixed count periods (1, 10, 20, 30, 40, 50, 60 or 100 s)
- Display of dose rate ($\mu\text{Sv/h}$) related to Cs-137
- GM tube voltage can be adjusted between 400 and 700 V, while showing count rate
- Max, min and average view for mSv/hr and CPS

Description

The GM tube

The Geiger tube is provided with a cap to protect the sensible window. When working with gamma radiation, the cap may remain in place. Remove the cap when alpha or beta radiation is involved.

Time

The counting period 1 s is achieved by selecting **CPS** (counts per second). All other periods, including infinity, is obtained with the **Pulse** setting. When the instrument has been turned off, the measurement period for **Pulse** is 10 s.

To change the period, press **HOLD ▲ HOLD** and set the measurement period with **▲** and **▼**.

The unit is seconds. The symbol "----" means infinite counting. Select with **CAL** / cancel with **HOLD**.

During **Pulse** measurements, the progress of each measurement period can be followed on the % -scale at the top.

GM voltage

The tube voltage can be adjusted in increments of 10 V. Nominal voltage is 500 V.

It is possible to follow the count rate (**cps**) in the secondary display, while the voltage is adjusted. Thus one can fairly easily estimate the knee voltage of the tube. This should be done by *increasing* the voltage as the response time is somewhat longer when the voltage is reduced.

Press **CAL** and adjust the voltage with **▲** and **▼**. Select with **CAL** / cancel with **HOLD**.

This setting returns to default value when the instrument is turned off.

Setting the default value for the GM voltage

Switch off with **OFF**. Keep **CAL** pressed down while turning on again (until the symbol **CAL** is displayed).

Set the voltage with **▲** and **▼**. Select with **CAL** / cancel with **HOLD**.

Statistics

The instrument keeps track of the maximum, minimum and average value of the **CPS** and **µSv/h**. The statistics functions start when the instrument is switched on.

The values are displayed like this:

HOLD ▼	Display: [HOLD	MAX]	Shows max. µSv/h and CPS	HOLD = exit
HOLD ▼▼	Display: [HOLD	MIN]	Shows min. µSv/h and CPS	HOLD = exit
HOLD ▼▼▼	Display: [HOLD	AVG]	Shows avg. µSv/h and CPS	HOLD = exit

Resetting

Instead of pressing **HOLD** in the end, press **CAL** to restart the statistics functions and use the current values for **CPS** and **µSv/h** as the new starting point.

Connections

Under a cover on the left side a connector that *looks like* a USB connector is found, which however *may not* be connected to USB devices. It is used for communication with external devices via special cables. It is possible to connect the instrument to the serial port on a PC – interface cable IM-159510. The same cable IM-159510 is used when connecting to a remote display IM 159810.

The built-in data logger

The instrument contains a data logger for 250 data sets. This enables measurements made in the field to be further processed on a computer. When the instrument is turned off, data are kept in the log – even during battery change.

Each data set contains these four values: **time**, **µSv/h**, **Mean CPS** and **Pulse**.

Data stored as the **Pulse** value depends on the logging mode and time interval for **Pulse**.

Manual logging: If **Pulse** is counting in finite time intervals, the number from the secondary display (last completed count) is stored.

If on the other hand **Pulse** is set up for infinite counting, the number from the primary display (ongoing count) is stored.

Automatic logging: **Pulse** data will accumulate – corresponding to infinite counting.

The hourglass icon turns on briefly to signal the logging of a data set.

If the log runs full, it will be indicated with a filled hourglass.

Manual logging

Quickly press **HOLD** twice to register a data set in the log.

Automatic logging

Press **HOLD ▲▲▲** (*display now shows free memory*) **HOLD**.

The display shows at the bottom [HOLD SET].

The primary display shows hours and minutes as "- :-" and the secondary display shows the seconds as ":-". This indicates that the automatic logging is *not* active.

With ▲ and ▼ you can now set the time interval to the following values:

1, 2, 5, 10, 20, 30, 40, 50 seconds; 1, 2, 5, 10, 20, 30, 40, 50 minutes; 1, 2, 3, 4, 5, 6 hours.

Select with **CAL** / cancel with **HOLD**. Pressing **CAL** erases any old data and starts logging.

When automatic logging is active, an empty hourglass icon is shown.

Resetting the log

If you want to erase the log – e.g. to make room for manual logging – carry out the above procedure, and set time as "- :- / :-".

Communication with PC

Communication requires a special cable (IM-159510). The cable connects to the PC serial port - eventually via a serial to USB adapter.

Communications protocol

Use a communication program like Hyper Terminal.

Use these parameters: 9600 baud, 8 data bits, parity "none", 1 stop bit, Xon / Xoff protocol.

Most commands to the unit consist of one letter, eventually followed by a number.

The possible responses from the device depend on the command.

V Reads the identification string from the instrument.

B Turns "GM Pulse" sound off or on.

D Reads the current measurement data: CPS, µSv/h, Mean CPS and Pulse – comma separated.

R Reads the current view from the primary display.

F Reads the operating mode of the instrument according to the table:

Returned value	1	2	3	4
Mode	CPS	µSv/h	Mean CPS	Pulse

G Reads the contents of the data logger in tabular form (comma separated), and stops any ongoing logging.

#07 Control code 07 (^G, ASCII "Bell") corresponds to a manual log command (two times HOLD)

(The command only works when the instrument is in REMOTE mode).

Optional accessories: Remote display and Cables

Remote display for Impo hand held instruments: IM-159810

Interface cable for PC or remote display: IM-159510

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© A/S Søren Frederiksen, Ølgod.