### Brugsvejledning til GM-sensor 5135.70

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2014.12.02

5135.70 AC

This GM Sensor is used for measuring radioactivity. The sensor detects alpha, beta and gamma radiation.

The sensor can be used with the following units: Electronic Counter 2002.50, Geiger-Müller Counter 5136.00 and Coincidence Box 5138.00

Furthermore the unit can be connected to any unit which can record TTL pulses and which can supply the sensor with a bias voltage of 5 VDC. See additional information under the Technical Specifications below.

The sensor is supplied with a mounting rod fitting the Bench for experiments with radioactivity 5141.00.

The sensor comprises a 500 V power supply and a GM tube mounted in the same housing.

The sensor is supplied with a power indicator LED. A sound generator in the sensor gives a "beep" whenever a radioactive particle is detected. This sound can not be switched off.


### IMPORTANT

The GM sensor is supplied with a window made of very thin and fragile material. The sensor is therefore supplied with a protective cover. The cover has a small air hole which must not be covered while the cover is put on and taken off the sensor. The cover must be removed when making measurements of alpha and beta radiation, while it can remain in place during measurements of gamma radiation.

### Technical Specifications:

The GM sensor can detect alpha, beta and gamma radiation.

Connecting jack:	6.3 mm jack
Supply voltage:	5 volts DC
Current drain:	ca. 10 mA
Output signal:	TTL
Size, L x B x H:	130 x 25 x 35 mm

### Window:

Thickness:	1.5 to 2.0 mg/cm <sup>2</sup>
Effective diam.:	9.14 mm
Material:	Mica

### Gas:

Neon and Halogen

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© A/S Søren Frederiksen, Ølgod