11.08.14

5141.00 AC

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
Bench for experiments with radioactivity

11.08.14

5141.00 AC

With this bench for radioactive sources and GM tubes a stable and suitable setup is achieved which allows you to perform reproducible experiments with radiation.

Three different types of sources may be used. See below for the relevant part numbers.

This equipment comprises

- Rail for mounting bench 40 cm long with mm scale
- Holder for radioactive sources with \varnothing 6 mm steel rod for hanging the absorber plates
- Rider for GM tube
- Box of absorber plates:
Lead plates (coated):
10 pcs. 1 mm, 5 pcs. 2 mm
Aluminium plates:
4 pcs. 0.5 mm, 6 pcs. 2 mm, 2 pcs. 3 mm
- Allen wrench for the adjustable distance indicator
- This instruction manual

Mounting absorber plates

First, place the GM tube and the source holder, with enough room for all the absorber plates you want to use at once. A distance of approx. 25 mm is appropriate for beta and gamma sources.

Next adjust the free length of steel rod to allow the absorber plates to be easily positioned without moving the tube or the source holder.

Alpha radiation has a very limited range in air. The distance in this case should therefore only be a few millimetres. Don't forget to remove the protective cap from the GM tube when working with an alpha source.

Distance Scale

The bench is fitted with a tape measure. Depending on the tube and source type, there will be a certain minimum distance between them, which is used to correct the raw distance values read from the scale. (It is recommended that you always note the raw distance readings, and perform the correction later.) The adjustable distance indicator might be moved using the supplied Allen wrench. Thereby one can obtain a value for the distance correction which is easier to use by mental arithmetic.

The following values apply only to Frederiksen's GM tube No. 5125.15 and Risø Sources No. 5100.00:

The distance from the GM tube's window to the front edge of the plastic housing is 4.5 mm.

The radioactive compound in Risø sources are located 3.0 mm behind the source tip.

With the source screwed tightly into the holder, the very source is pulled another 1.5 mm back (approx.). This means that when the GM tube is very close to (touching) the source holder, there is a distance of approx. 9.0 mm between the tube window and the source's centre.

Part number for bench

5141.00

5141.20

5141.10

Corresponding source type

Risø / Frederiksen

25 mm disc

12 mm cylindrical

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© A/S Søren Frederiksen, Ølgod